

Model: DMS 300-4L

Multi-Parameter: ECG, Respiration, ABP, and Body Position



- ECG Memory transfer from either USB cable interface or mini SanDisk card.
- Memory \geq 4GB
- Record multi-day Holter ECG and overnight Respiration.
- Records patient body position.
- HDMI for artifact isolation.
- 3- and 12-lead ECG recordings (simply plug-in the appropriate ECG cable).
- Cardiac Event trans-telephonic ECG for sending of symptomatic ECG's.
- Separate SAECG file recording at 1,024 Hz and 16-bit.
- LCD: Real-time ECG at any time, Heart Rate, Time, Battery status, ECG cable status, Respiration rate.
- Independent pacemaker spike channel with detection of mini pacemaker spikes and 10,240 sample rate.
- 16 bits A/D resolution with frequency response of 0.05 to 150 Hz.
- 5 day Holter recordings.

- 1 AAA battery power supply.
- High speed USB interface.
- Sample rate: read 4,096 Hz, write 128 Hz.
- Implantable Cardioverter-Defibrillator (ICD) **Protection**.
- Digital reliability with 2 year warranty.

